# Immuglo™ Anti-Skin Antibody (IC/BMZ) Test System Monkey/Guinea Pig Esophagus Sections

For *in vitro* Diagnostic Use **IVD**

CLIA Complexity: High  
CDC Analyte Identification Code: 0449  
CDC Test System Identification Code: 28313

**REF** Code: 1104 48 determinations

## PRODUCT INSERT

### INTENDED USE

An indirect immunofluorescence antibody test for the detection and quantitation of anti-skin (anti-intercellular and anti-basement membrane) antibodies in human serum.

### SUMMARY AND EXPLANATION

The detection of anti-skin antibodies aids in the diagnosis, and sometimes prognosis, of chronic vesicular-bullous diseases including pemphigus, pemphigoid, cicatricial pemphigoid, and epidermolysis bullosa acquisita (EBA)<sup>1</sup>. Epithelial intercellular antibodies are diagnostic for pemphigus and occur in over 90% of sera from all active forms<sup>2</sup>. Because of the species specificity of IC antibodies, the use of the dual monkey/guinea pig esophagus substrates further distinguishes the intercellular antibodies of pemphigus vulgaris and vegetans, which preferentially react with monkey esophagus sections, from the antibodies of pemphigus foliaceus and erythematosus which react more strongly with guinea pig esophagus sections<sup>3-6</sup>. Antibodies to basement membrane antigens of stratified squamous epithelium occur in about 70% of active bullous pemphigoid, 50% of vesicular pemphigoid and EBA and 10% of cicatricial pemphigoid patients. The basement membrane antibodies usually react with both monkey esophagus and guinea pig esophagus sections<sup>3-6</sup>.

### PRINCIPLES OF PROCEDURE

In the indirect immunofluorescence method used in this kit, patients' sera are incubated on monkey/guinea pig esophagus sections to allow binding of antibodies to the substrate. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgG class are detected by incubation of these sections with fluorescein-labeled, anti-human IgG. Reactions are observed under a fluorescence microscope equipped with appropriate filters.

The presence of anti-intercellular and anti-basement membrane antibodies is demonstrated by an apple green fluorescence of these respective histologic structures. The titer (the reciprocal of the highest dilution giving a positive reaction) is then determined by testing serial dilutions<sup>2</sup>.

### PRODUCT INFORMATION

#### Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

#### Materials provided

**REF** Code: 1104 48 determinations

8 x	<b>SORB SLD 6</b>	6 well Substrate Slides, Monkey/Guinea Pig Esophagus
1 x 0.5 ml	<b>CONTROL + IC</b> *	Anti-IC Positive Control. Contains human serum.
1 x 0.5 ml	<b>CONTROL + BMZ</b> *	Anti-BMZ Positive Control. Contains human serum.

1 x 0.5 ml	<b>CONTROL -</b> *	Negative Control. Contains human serum.
1 x 5 ml	<b>CONJ FITC</b> *	Anti-human FITC Conjugate. Protect from light.
1 x 5 ml	<b>CONJ FITC EB</b> *†	Anti-human FITC Conjugate containing Evan's Blue. Protect from light.
1 x 60 ml	<b>BUF</b> *	Buffered Diluent.
2 vials	<b>BUF WASH</b>	Phosphate Buffered Saline (PBS). Dissolve each vial to 1 liter.
1 x 5.0 ml	<b>MOUNTING MEDIUM</b> *	Mounting Medium. Do not freeze.
1 x 1.0 ml	<b>EVANS</b>	Evan's Blue Counterstain.
1 x 12	<b>COVER SLD</b>	Coverslips.
<b>Optional Components</b>		
1 x 5 ml	<b>IgG-Conj FITC PA</b> *†	Anti-human FITC Conjugate containing Evan's Blue. Protect from light.

\* Contains < 0.1% NaN<sub>3</sub>

† Replaces conjugate without counterstain in Code numbers containing "EB"

#### Material required but not provided

Fluorescence microscope  
 Micropipette or Pasteur pipette  
 Serological pipettes  
 Staining dish (e.g. Coplin jar)  
 Small test tubes (e.g. 13 x 75 mm) and test tube rack  
 Distilled or deionized water  
 1 liter container  
 Wash bottle  
 Paper towels  
 Incubation chamber

#### WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>7</sup>.

WARNING - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

#### SPECIMEN COLLECTION AND PREPARATION

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

#### PROCEDURE

##### Test Method

#### A. Screening

1. Dilute each patient serum 1:10 with the Buffered Diluent provided (10  $\mu$ l serum + 90  $\mu$ l Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50  $\mu$ l of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50  $\mu$ l) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50  $\mu$ l) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evans blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

#### B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following **Steps 5 through 13** to determine the titer. Each test run should include the one slide with the undiluted Negative Control and undiluted, 1:2, 1:4, 1:8 and 1:16 dilutions of the Anti-Intercellular Positive Control. Make serial two-fold dilutions of the patient's serum starting at 1:10 (see following diagram). Using another slide, a serum positive for intercellular antibodies may be tested at dilutions ranging from 1:10 to 1:320. An additional slide should be used to obtain endpoints for those patients' sera still positive at a 1:320 dilution.

A serum positive for anti-basement membrane antibodies when screened at doubling dilutions of 1:10 through 1:80 may be reported as positive,  $\geq 80$ . If instead, a serum has been screened at 1:20 and 1:40, then prepare additional two fold dilutions of 1:80 to 1:640 and use an additional slide. Include the anti-basement membrane control in one well and the Negative Control in another well. **DO NOT** dilute these controls.

The reciprocal of the highest dilution producing a positive reaction is the titer.

### Preparation of Serial Dilutions

For sera positive for intercellular antibodies upon screening, number twelve tubes 1 through 12. Add 0.9 ml of buffered diluent to tube 1 and 0.2 ml to tubes 2 through 12. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4	5	6, etc.
Serum	0.1 ml					
	+					
Buffered Diluent	0.9 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
		↻	↻	↻	↻	↻
Transfer		0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
Final dilution	1:10	1:20	1:40	1:80	1:160	1:320, etc.

For sera positive for basement membrane antibodies when tested at 1:20 and 1:40 dilutions, number seven tubes 1 through 7 and make dilutions as described above.

### Quality Control

Both the Positive and Negative Controls should be included with each test run. The Negative Control should show no specific fluorescence of the intercellular or basement membrane areas. The Anti-Intercellular Positive Control and the Anti-Basement Membrane Positive Control both should have 2+ or greater staining intensity of the epithelial intercellular and the basement membrane areas, respectively, on both substrates. Upon titration of the anti-intercellular positive control, an endpoint titer of  $8 \pm$  one doubling dilution should be obtained. Sample reactions are provided in Figures 1 and 2 at the end of this document.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control.
- Problems with the optical system of the fluorescence microscope. These may include: improper alignment, use of the bulb beyond the expected performance life, etc.
- Allowing the slide to dry during the procedure.
- Improper preparation of serial dilutions of controls

### RESULTS

The results of the tests for anti-intercellular antibodies should be reported as negative (<10) or positive with titer. The results of the tests for anti-basement membrane antibodies should be reported as negative (<10), positive ( $\geq 80$  or  $\geq 640$ ) or alternatively, positive with titer.

Read for specific staining of the epithelial intercellular substance or basement membrane zone. Various other tissue antibodies such as ANA, mitochondrial, smooth muscle, endomysial, basal cell cytoplasmic and upper cytoplasmic antibodies and antibodies to striated muscle may also be observed on monkey and guinea pig esophagus sections<sup>4</sup>. Sera giving any of these reactions should be reported as negative for antibodies to intercellular and to basement membrane antigens. Any sera giving nuclear staining reactions may be tested on Antinuclear Antibody (ANA) Test (HEp-2 Cells) *Catalog Nos. 1102 or 1103*, Antinuclear Antibody (ANA) Test (Mouse Liver Sections) *Catalog Nos. 1100 or 1101*. Any sera giving smooth muscle or mitochondrial staining reactions may be tested on Autoantibody Test System (Mouse Kidney/Stomach Sections) *Catalog No. 1107*.

#### LIMITATIONS OF THE PROCEDURE

In some sera both intercellular and basement membrane antibodies may be present in varying titers as indicated by different staining patterns. Both antibody patterns and their titers should be reported. Occasionally sera may exhibit strong staining for ANA, antibodies to the cytoplasm, and non-specific staining of spinous or basal cell layer of epithelium. These may interfere with the ability to detect anti-intercellular and anti-basement membrane antibodies, respectively. In such cases titrating the serum may permit visualization of the skin antibodies. In other cases they may be in lower titer than the ANA or cytoplasmic antibodies.

In the latter situation, the serum should be reported as negative for intercellular and basement membrane antibodies, but it should be mentioned that skin antibodies cannot be excluded. In cases such as these, biopsy studies by direct immunofluorescence are recommended.

The goat anti-human IgG FITC Conjugate supplied in this kit is primarily heavy chain specific but has some light chain activity. It reacts primarily with IgG class autoantibodies, but may, to a lesser degree, react with light chains of other classes such as IgM.

A positive basement membrane antibody reaction with a linear pattern is consistent with a diagnosis of pemphigoid (i.e. bullous pemphigoid or other forms of pemphigoid) or of epidermolysis bullosa acquisita. The indirect immunofluorescence method using these substrates cannot distinguish between these two diseases.

A positive intercellular reaction is not diagnostic of pemphigus. Such reactions may also be due to blood group antibodies or pemphigus-like antibodies. It is recommended that all sera positive for intercellular antibodies be absorbed to distinguish those sera giving positive intercellular staining due to blood group antibodies<sup>3</sup>.

Therefore, clinicians need to consider the serum findings together with the clinical findings and other laboratory studies, notably histopathology and direct immunofluorescence of skin biopsy specimens, towards the diagnosis of pemphigus or pemphigoid.

#### EXPECTED VALUES

Almost all pemphigus patients have IC antibodies. The use of the dual substrates increases the sensitivity of detecting intercellular antibodies as indicated by the data in Table 1 at the end of this document on 123 pemphigus patients.

In addition to increasing the sensitivity, the use of both tissue substrates enables the immunologic differentiation of pemphigus vulgaris from pemphigus foliaceus in most cases<sup>4-6</sup>. As illustrated in Table 2, pemphigus vulgaris sera generally give higher titers and/or brighter staining on monkey esophagus than on guinea pig esophagus, while the opposite is true for pemphigus foliaceus sera.

Comparison of titers of pemphigus antibodies in previous and current sera from patients with a proven diagnosis of pemphigus affords valuable prognostic information since the fluctuations in titer tend to parallel the disease activity as shown in Figure 4. Titers afford useful information for the adjustment of doses of therapeutic agents; patients with residual titers but no lesions may be expected to have relapses if therapy is withdrawn<sup>6</sup>.

The percentage of pemphigoid patients with anti-basement membrane antibodies varies with the form of the disease as shown by Table 3.

Differentiation of the antibodies of epidermolysis bullosa acquisita and pemphigoid requires special test procedures<sup>9-10</sup>. The recommended method is an indirect immunofluorescence test of basement membrane antibodies on normal skin split with 1M NaCl<sup>9</sup>. Titers of basement membrane antibodies usually do not parallel disease activity as do pemphigus antibodies.

#### PERFORMANCE CHARACTERISTICS

Sera obtained from 38 normal subjects, 10 confirmed cases of pemphigus vulgaris, 10 confirmed cases of pemphigus foliaceus and 20 confirmed cases of pemphigoid or epidermolysis bullosa acquisita were tested on this kit and another commercially available skin antibody kit. Sera were tested according to the procedures recommended by the manufacturers. Results using the Immuglo™ Anti-Skin Antibody (IC/BMZ) Test System showed increased sensitivity as well as comparable specificity. See Table 4.

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
#### REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI


1. Beutner EH, Chorzelski TP and Kumar V (Eds) "Immunopathology of the Skin", 3rd Ed, John Wiley and Sons, New York, 1987.
2. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
3. Kumar, V and Beutner EH, Monkey esophagus: A unique antigenic substrate in immunodermatology. In "immunopathology of the skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 65-89, 1987.
4. Sabolinski, ML, Beutner EH, Krasny S, Kumar V and Chorzelski TP. Substrate specificity of anti-epithelial antibodies of pemphigus vulgaris and pemphigus foliaceus sera in immunofluorescence tests on monkey and guinea pig esophagus sections. *J Invest Dermatol* 88: 545-549, 1987.
5. Beutner EH, Sabolinski ML, Huang J, Krasny SA, Kumar V and Chorzelski TP. Indirect immunofluorescent differentiation of autoantibodies of pemphigus vulgaris and pemphigus foliaceus. *J Invest Dermatol* 86: 463, 1986.
6. Krasny S, Beutner EH and Chorzelski TP. Specificity and sensitivity of indirect and direct immunofluorescence findings in the diagnosis of pemphigus. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 207-247, 1987.
7. Biosafety in Microbiological and Biomedical Laboratories Centers for Disease Control, National Institutes of Health, (HHS Pub. No{CDC} 93-8395) 1999.
8. Gammon WR, Wilson BD and Briggaman RA. Epidermolysis bullosa acquisita. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 383-398, 1987.
9. Medicamojslovic L, Fenske W and Espinoza C. Epidermolysis bullosa acquisita - direct immunofluorescence and ultrastructural studies. *Am J Dermatopathol* 9: 324-333, 1987.
10. Gammon WR, Kowalewski C, Chorzelski TP, Kumar V, Briggaman RA and Beutner EH. Direct immunofluorescence on sodium chloride-separated skin in the differential diagnosis of bullous pemphigoid and epidermolysis bullosa acquisita. *J Am Acad Dermatology*, 1989.

**IC staining for IgG in Pemphigus**

**Monkey Esophagus**

Figure 1

**Guinea Pig Esophagus**

Figure 2

**BMZ staining for IgG in Pemphigoid**


Figure 3


**Table 1: Sensitivity of Indirect Immunofluorescence Test for Antibodies On Different Substrates**

Tested on	Negative	Positive	% Positive
Monkey Esophagus alone	12	111	89
Guinea Pig Esophagus alone	20	103	81
Together	1	122	99

**Table 2: Differentiation of Pemphigus Vulgaris from Pemphigus Foliaceus by Indirect Immunofluorescence Test**

Higher Titer and/or Brighter Staining on:	P. Vulgaris Sera	P. Foliaceus Sera
Monkey Esophagus	73	0
Guinea Pig Esophagus	0	25
No Difference	2	5

**Figure 4: Pemphigus Antibody Titers and Disease Activity**  
 Reprinted with permission from: Krasny S, Beutner EH and Chorzelski TP<sup>6</sup>

**Table 3: Incidence of Basement Membrane Antibodies Detected by Indirect Immunofluorescence on Monkey Esophagus Substrate**

Clinical Condition	% Positive
Bullous Pemphigoid	70
Cicatricial Pemphigoid and Brunsting-Perry Variant	10
Vesicular Pemphigoid	50
Epidermolysis Bullosa Acquisita	50
Normal Controls	0

**Table 4: Comparisons of Kits Using Monkey/Guinea Pig Esophagus vs. Monkey Esophagus Substrate Alone for the Detection of Anti-IC and Anti-BMZ Antibodies by Indirect Immunofluorescence**

Clinical Condition	No. of Sera	% Positive		
		IMMCO Monkey Esophagus	IMMCO Guinea Pig Esophagus	OTHER Monkey Esophagus
Pemphigus Vulgaris	10	10/100	5/50	2/20
Pemphigus Foliaceus	10	6/60	10/100	3/30
Pemphigoid/ EBA	20	20/100	ND	17/85
Normal Controls	38	4*/11	0	0

EBA = Epidermolysis bullosa acquisita  
 ND = Not Done

\*These four sera upon treatment with absorbent were negative for intercellular staining on monkey esophagus sections due to blood group antibodies.

*For technical assistance please contact:*

**IMMCO**  
DIAGNOSTICS

**IMMCO Diagnostics, Inc.**  
60 Pineview Drive  
Buffalo, NY 14228-2120  
Telephone: (716) 691-0091  
Fax: (716) 691-0466  
Toll Free USA/Canada: 1-800-537-TEST  
E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé  
EMERGO Group, Inc.  
Molenstraat 15, 2513 BH, The Hague,  
The Netherlands  
Tel (+31) 345 8570, Fax (+31) 346 7299  
[www.emergogroup.com](http://www.emergogroup.com)

REV.APR2009

Document No. PI4108CE