ImmuLisa™ Anti-Cardiolipin Antibody (ACA) Screen ELISA

IVD For in vitro diagnostics use

CLIA Complexity: High
 CDC Analyte Identification Code: 0434
 CDC Test System Identification Code: 28515

PRODUCT INSERT

REF 1118S Anti-Cardiolipin Antibody Screen ELISA 96 Determinations

INTENDED USE

An enzyme linked immunoassay (ELISA) for the detection of IgG, IgA and IgM class anticardiolipin antibodies as an aid in assessing the risk of thrombosis in individuals with anti-phospholipid syndrome.

SUMMARY AND EXPLANATION

Antiphospholipid antibodies are a heterogeneous group of autoantibodies against negatively charged phospholipids¹. They are detected primarily by the anti-cardiolipin antibody (ACA) test, the biological false positive test for syphilis and the lupus anticoagulant test. These three tests detect related, but not necessarily identical antibodies. Thus, more than one of these tests is sometimes necessary to identify antiphospholipid antibodies.

Of the various tests for the detection of antiphospholipid antibodies, the anti-cardiolipin antibody test performed by ELISA is the most sensitive². The presence of anti-cardiolipin antibodies helps to identify patients at risk of venous and/or arterial thrombosis often accompanied by thrombocytopenia, a syndrome referred to as antiphospholipid syndrome¹⁻¹². The antiphospholipid syndrome most commonly occurs in patients with systemic lupus erythematosus (SLE) or lupus-like disease where the criteria for SLE are not fulfilled⁵⁻⁷. High levels of anti-cardiolipin antibodies occur in thrombosis, fetal loss, thrombocytopenia and several other disorders¹⁻¹⁵. Low levels of anti-cardiolipin antibodies are found in a variety of clinical disorders which are unrelated to antiphospholipid syndrome. Therefore, low levels of these antibodies are of limited significance.

IgG and IgA class anti-cardiolipin antibodies appear to be more closely associated with antiphospholipid syndrome than the IgM class antibodies. However, IgM antibodies appear to be more influenced by treatment^{5,10}. Low levels of IgM antibodies can be identified in other autoimmune diseases such as rheumatoid arthritis, primary Sjögren's Syndrome, drug induced lupus erythematosus, Lyme disease, and syphilis^{8,10}.

The ImmuLisa™ ACA Screen ELISA offers a convenient method of screening patient sera for the presence of ACA before further identifying positive samples for individual isotypes with the ImmuLisa™ ACA-IgG, -IgA and-IgM ELISA.

PRINCIPALS OF THE PROCEDURE

Cardiolipin antigen is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrator and diluted patient sera are added to separate wells, allowing any cardiolipin antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human conjugate is added to each well. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the appropriate class. After washing away any un-bound conjugate, specific enzyme substrate (pNPP) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in ELISA units per milliliter (EU/ml).

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials¹⁶. **WARNING** - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this kit insert to en-sure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc.

Use good laboratory techniques to minimize microbial and chemical contamination. Do not use after expiration date.

Materials provided

ImmLISA™ Anti-Cardiolipin Antibody Screen ELISA **REF** 1118S

Kits contain sufficient reagents to perform 96 determinations each.

| | | |
|------------|--|--|
| 12 x 8 | MICROPLATE ACA | Microplate with individual breakaway microwells coated with cardiolipin antigen. |
| 1 x 1.5 ml | CALIBRATOR ACA SCREEN * | Ready to use Calibrator (<i>green cap</i>). Human serum containing antibodies to cardiolipin. |
| 1 x 1.5 ml | CONTROL + ACA SCREEN * | Ready to use Positive Control (<i>red cap</i>). Contains human serum positive for ACA. |
| 1 x 1.5 ml | CONTROL - * | Ready to use Negative Control (<i>white cap</i>). Contains human serum. |
| 1 x 12 ml | IgA/IgG/IgM - CONJ ALKPHOS * | Ready to use anti-human Alk. Phos. Conjugate . Color coded pink. |
| 2 x 60 ml | DIL * | Ready to use Serum Diluent . Color coded blue. |
| 1 x 12 ml | SUBSTRATE * | Ready to use Enzyme Substrate . Contains pNPP. Protect from light. |
| 1 x 12 ml | STOP | Ready to use Stop Solution . |
| 2 x | BUF WASH | Powder Wash Buffer . Reconstitute to one liter each. |

* Contains < 0.1% NaN₃

Symbols used on labels:

LOT Lot number

REF Catalog number

 Use by

 Storage temperature

 Read instructions for use

IVD In vitro diagnostic use

 Manufacturer

 Number of Tests

EN

Materials Required But Not Provided

- Pipette capable of dispensing 5 μ l and 1000 μ l
- Disposable pipette tips
 - Clean test tubes 12 x 75 mm and test tube rack
 - Deionized or distilled water
 - Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
 - Squeeze bottle to hold diluted wash buffer
 - Timer
 - Absorbent paper
- Automated microplate washer capable of dispensing 200 μ l

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

PROCEDURE

Procedural Notes

- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-26°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

Assay procedure

1. **ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-26°C) PRIOR TO BEGINNING THE ASSAY.**
2. Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.
3. Prepare a **1:201** dilution of the patient specimen by mixing **5 μ l** of the patient specimen with **1.0 ml** of Serum Diluent.
4. Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator.

QUALITATIVE DETERMINATION

5. Add **100 µl** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record.
Note: Include one well with **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of this well should not be greater than 0.3.
6. Incubate **30 minutes** (± 5 minutes) at room temperature on a level surface.
7. Wash step: Thoroughly aspirate the contents of each well. Add 200-300µL of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
8. Add 1 00µL of the Conjugate to each well.
9. Incubate the wells for **30 minutes** (± 5 min) at room temperature.
10. Wash step: Repeat step 7.
11. Add 100µL of En yme Substrate to each well.
12. Incubate for **30 minutes** (± 5 min) at room temperature.
13. Add 1 00µL of Stop Solution to each well. Maintain the same se uence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
14. Read the absorbance (OD) of each well at 405nm using a single or dual wavelength microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrator, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <20 EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml.

RESULTS

Calculations

The following method should be used to determine whether the specimen is positive or negative for ACA.

QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator}} \times \text{EU/ml of Calibrator} = \text{EU/ml Test Sample}$$

Interpretation

The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

| ACA Screen Value | Interpretation |
|------------------|----------------|
| <20 EU/ml | Neg (-) |
| 20-25 EU/ml | Borderline |
| >25 EU/ml | Pos (+) |

The literature suggests that, low positive anti-cardiolipin antibody levels may occur in a variety of clinical disorders unrelated to antiphospholipid antibody syndrome. Hence according to the investigators recommendations the diagnosis of antiphospholipid antibody syndrome should be made only when the test results are moderately or highly positive¹⁴.

LIMITATIONS OF THE PROCEDURE

A diagnosis cannot be made on the basis of ACA results alone. The results of other laboratory tests and clinical findings must be considered. Tests for rheumatoid factor (RF) should be performed because RF may interfere with this assay. When a negative ACA test occurs in the presence of clinical indications, a lupus anticoagulant test or other additional testing is indicated. ACA also occur transiently in a variety of infectious diseases. In these cases patients positive for ACA should be retested following an appropriate interval. Confirmed active or seropositive syphilis patients can have elevated ACA levels. To rule out syphilis, confirmatory tests should be performed.

EXPECTED VALUES

The incidence of ACA in various disease conditions is summarized in the following tables:

Incidence of ACA in SLE^{15,17}

| Antibody Isotype | % Incidence |
|------------------|-------------|
| IgG | 39-44 |
| IgA | 17-57 |
| IgM | 5-33 |
| Any isotype | 53 |

Disease Association of ACA

| Condition | % Incidence |
|-----------------------------|-------------|
| Recurrent Venous Thrombosis | 28-71 |
| Recurrent Fetal Loss | 28-64 |
| Transverse Myelitis | 50 |
| Hemolytic Anemia | 38 |
| Thrombocytopenia | 27-33 |
| Arterial Occlusions | 25-31 |
| Livedo Reticularis | 25 |
| Pulmonary Hypertension | 20-40 |

PERFORMANCE CHARACTERISTICS**Comparison of ELISA Methods for the Detection of ACA**

| | | Immco™ IgG-ACA | | |
|-------------|---------|----------------|---------|---------|
| | | POS (+) | NEG (-) | TOT (=) |
| ImmuLisa™ | POS (+) | 31 | 0 | 31 |
| ACA IgG/A/M | NEG (-) | 2 | 42 | 44 |
| TOT (=) | | 33 | 42 | 75 |

relative specificity: 95.5%
 relative sensitivity: 100.0%
 relative agreement: 97.3%

Precision:

Two samples with known concentrations of ACA were assayed in 10 replicates over a period of two weeks. Intra- and inter-assay coefficient of variation (CV) were as follows:

| | Intra- plate CV | Inter-plate CV |
|----|--------------------|-------------------|
| 1. | 7.6% | 5.9% |
| 2. | 10.4% | 6.3% |

Recovery:

Samples with known concentrations of ACA were mixed with the appropriate dilutions of another positive sample with known amounts of ACA. The ACA values of the spiked samples were then determined and from the values obtained the percent recovery calculated.

| | ACA added | ACA obtained | % Recovery |
|----|-----------|--------------|------------|
| 1. | 127.0 | 135.8 | 106.9 |
| 2. | 103.8 | 103.2 | 99.4 |
| 3. | 64.3 | 61.2 | 95.2 |

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REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPÜIE • BIBLIOGRAFIA

1. Harris EN, Phil M, Asherson RA et al. Antiphospholipid antibodies - antibodies with a difference. *Ann Rev Med*; 1988, 39:261-271.
2. Harris EN, Gharavi AE, Wasley GD and Hughes GRV. Use of an enzyme-linked immunosorbent assay and inhibition studies to distinguish between antibodies to cardiolipin from patients with syphilis or autoimmune disorders. *J Infect Dis*; 1988, 157:23-31.
3. Lockshin MD. Anti-cardiolipin antibody. *Arth Rheum*; 1987, 30:471-472.
4. Derksen RHW, Beisma D, Bouma BN et al. Discordant effects of prednisone on anticardiolipin antibodies and the lupus anti-coagulant. *Arth Rheum*; 1986, 29:1295-1296.
5. Alarcón-Segovia D, Deleze M, Oria CV et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. *Medicine*; 1989, 68:353-365.
6. Triplett DA. Antiphospholipid antibodies and recurrent pregnancy loss. *Am J Reprod Immunol*; 1989, 20:52-67.
7. Montecucco C, Longhi M, Caporali R et al. Hematological abnormalities associated with anticardiolipin antibodies. *Hematologica*; 1989, 74:195-204.
8. Asherson RA, Khamashta MA, ORDI-Ros J et al. The "primary" antiphospholipid syndrome: major clinical and serological features. *Medicine*; 1989, 68:366-374.
9. Walport MJ. Pregnancy and antibodies to phospholipids. *Ann Rheum Dis*; 1989, 48:795-797.
10. Machworth-Young CG, Loizou S and Walport MJ. Antiphospholipid antibodies and disease. *Medicine*; 1989, 72:767-777.
11. Harris EN, Gharavi AE and Hughes GRV. Antiphospholipid antibodies. *Clin Rheumatol*; 1985, 1:591-609.
12. Lubbe WF, Palmer SJ, Butler WS et al. Fetal survival after prednisone suppression of maternal lupus anticoagulant. *Lancet*; 1983, 1:1361-1363.
13. Hogan MJ, Brunet DG, Ford PM et al. Lupus anticoagulant, antiphospholipid antibodies and migraine. *Can J Neurological Sci*; 1988, 15:420-425.
14. Harris EN and Pierangeli S. Antiphospholipid antibodies: method of detection. *Am J Reproduct Immunol*; 1992, 28:208-210.
15. Silveira LH, Lopez LR, Uzzle RT et al. IgA, IgG and IgM anticardiolipin antibodies in black patients with systemic lupus erythematosus. *Arthr Rheum 56th Annual Meeting*; 1992, S125.
16. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health; 1993, (HHS Pub. No [CDC] 93-8395).
17. Alarcón-Segovia D, Delezé M, Oria CV et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. A perspective analysis of 500 consecutive patients. *Medicine*; 1987, 68:353-365.

For technical assistance please contact:

IMMCO Diagnostics, Inc.

60 Pineview Drive

Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immco.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,
The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com