IMMCO
DIAGNOSTICS

ANA/AMA/ASMA/AGPA/LKM Antibody Test System

IVD

PRODUCT INSERT

REF 1134LKM HEp-2/Mouse Liver/Kidney/Stomach 48 Determinations

REF 1136C Mouse Liver/Kidney/Stomach 48 Determinations

REF 2152-3 Mouse Liver/Kidney/Stomach Slide 8 Wells

REF 2190-LKM HEp-2/Mouse Liver/Kidney/Stomach Slide 8 Wells

INTENDED USE

Indirect immunofluorescence (IF) antibody tests for the detection and quantitation of antinuclear antibodies (ANA), anti-mitochondrial antibodies (AMA), anti-smooth muscle antibodies (ASMA), anti-gastric parietal cell antibodies (AGPA), and anti-liver/kidney/microsomal (LKM) antibodies in human serum.

SUMMARY AND EXPLANATION

Antinuclear antibodies (ANA) detected by indirect immunofluorescence, aid in the diagnosis of connective tissue disorders including systemic lupus erythematosus (SLE), mixed connective tissue disease, Sjögren's syndrome and scleroderma¹⁻⁵. ANA occur in about 95% of SLE patients as well as patients with other connective tissue diseases. ANA may also occur in other disorders such as chronic active hepatitis and primary biliary cirrhosis⁶⁻⁸.

Anti-mitochondrial antibodies (AMA) occur in over 90% of primary biliary cirrhosis cases, 3-11% of chronic active hepatitis patients and are absent in patients with extrahepatic biliary obstruction and in other liver diseases. The universal presence of anti-mitochondrial antibodies in primary biliary cirrhosis and their virtual absence in extrahepatic jaundice makes their detection of considerable value in the differential diagnosis⁶⁻¹².

Anti-smooth muscle antibodies (ASMA) in high titer (>160) occur in the majority of cases of chronic active hepatitis and in intermediate titers (40-80) in acute viral hepatitis. Occasionally they may occur in cases of primary biliary cirrhosis where they are also found in intermediate titers. The significance of titers of 20-40 is doubtful since these titers may occur in normal individuals^{13,14}.

Anti-gastric parietal cell antibodies (AGPA) are commonly associated with pernicious anemia and chronic atrophic gastritis where they occur in about 90% and 50% of cases, respectively. However, they are not disease specific as they may occur in low frequency in other disorders. Although healthy individuals may have gastric parietal cell antibodies, this finding may reflect asymptomatic atrophic gastritis. Negative findings for gastric parietal cell antibodies provide strong evidence for excluding pernicious anemia¹⁵⁻¹⁷.

Anti-liver/kidney/microsomal antibodies (LKM) are microsomal antibodies that exhibit characteristic reaction pattern staining of the cytoplasm of liver and proximal tubules of the kidney tissue. They are found in a group of patients with autoimmune hepatitis^{18,20}. LKM antibodies are primarily of the IgG isotype and can be differentiated by various reaction patterns.

PRINCIPLES OF PROCEDURE

In the indirect immunofluorescence method used in this kit, patients' sera are incubated on HEp-2/Mouse liver/kidney/stomach sections to allow binding of antibodies to the substrate. Unbound antibodies are removed by rinsing the slide. Bound antibodies of the IgG class are detected by incubation of the substrate with a fluorescein-labeled, anti-human IgG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of ANA, ASMA, AMA, AGPA and LKM antibodies are demonstrated by an apple green fluorescence of specific histologic structures in the substrate. The titers (the reciprocal of the highest dilution giving a positive reaction) are determined by testing serial dilutions²¹.

PRODUCT INFORMATION

Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

EN

Materials provided

ImmuGlo™ HEp-2/Mouse Liver/Kidney/Stomach Kit **REF** 1134LKM

ImmuGlo™ Mouse Liver/Kidney/Stomach Kit **REF** 1136C

Kit contains sufficient reagents to perform 48 determinations

6 x	SORB SLD 8	8 well Substrate Slides, HEp-2/Mouse Liver/Kidney/Stomach (1134LKM)
6 x	SORB SLD 8	8 well Substrate Slides, Mouse Liver/Kidney/Stomach (1136C)
1 x 0.5 ml	CONTROL + ANA *	ANA Positive Control. Contains human serum. (1134LKM)
1 x 0.5 ml	CONTROL + AMA *	AMA Positive Control. Contains human serum. (1134LKM, 1136C)
1 x 0.5 ml	CONTROL + LKM	LKM Positive Control. Contains human serum. (1136C)
1 x 0.5 ml	CONTROL - *	Negative Control. Contains human serum.
1 x 5 ml	IgG-CONJ FITC *	Anti-human IgG FITC Conjugate. Protect from light
1 x 5 ml	IgG-CONJ FITC EB *†	Anti-human IgG FITC Conjugate containing Evan's Blue. Protect from light.
1 x 60 ml	BUF *	Buffered Diluent.
2 vials	BUF WASH	Phosphate Buffered Saline (PBS). Dissolve each vial to 1 liter.
1 x 5.0 ml	MOUNTING MEDIUM *	Mounting Medium. Do not freeze.
1 x 1.0 ml	EVANS	Evan's Blue (EB) Counterstain.
1 x 12	COVER SLD LONG	Long Coverslips (1134LKM)
1 x 12	COVER SLD	Coverslips (1136C)

* Contains < 0.1% NaN₃

† Replaces conjugate without counterstain in Code numbers containing "EB"

Symbols used on labels:

LOT Lot number

REF Catalog number

 Use by

 Storage temperature

 Read instructions for use

IVD In vitro diagnostic use

 Manufacturer

 Number of Tests

Material required but not provided

- Fluorescence microscope

EN

- Micropipette or Pasteur pipette
- Serological pipettes Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation chamber

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials¹⁹. **WARNING** - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center. Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc.. Do not use beyond expiration date.

SPECIMEN COLLECTION AND PREPARATION

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

PROCEDURE

Test Method

A. Screening

1. Dilute each patient serum 1:10 with the Buffered Diluent provided (10 µl serum + 90 µl Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 µl) of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the

slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 μ l) to each well.

9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If optional conjugate without counterstain is used (see optional components in Materials Provided Section), 2-3 drops of Evan's Blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:10. The reciprocal of the highest dilution producing a positive reaction is the titer.

Preparation of Serial Dilutions

Number six tubes 1 through 6. Add 0.9 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4	5	6
Serum	0.1 ml					
	+					
Buffered Diluent	0.9 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
		⇌	⇌	⇌	⇌	⇌
Transfer		0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
Final dilution	1:10	1:20	1:40	1:80	1:160	1:320 etc.

QUALITY CONTROL

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence. The ANA Positive Control should have 2+ or greater staining intensity of the nuclei of the kidney with a predominantly homogeneous pattern. If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

INTERPRETATION OF RESULTS

The results of the tests for ANA, AMA, ASMA, AGPA, and anti-LKM antibodies should be reported as negative (with titer less than 10), positive (with titer greater than or equal to 320) or, alternatively, positive with specific endpoint titer. Read only fields which contain specific staining of the nuclei for ANA, the kidney tubules for AMA, the kidney blood vessel walls for ASMA, gastric parietal cells only for AGPA, and liver parenchymal cells and kidney proximal tubules for LKM antibodies. All other reactions should be reported as negative for ANA, AMA, ASMA, AGPA and/or LKM. ANA can be detected on all substrates but should be quantified on the kidney or HEp-2 cells. The nuclear staining patterns observable with the kidney substrate or HEp-2 cells provided include homogeneous, peripheral (rim), speckled and nucleolar. The centromere staining pattern (including mitotic figures) is seen most easily on HEp-2 cells. These nuclear staining patterns are described below. They may be one or a combination of several staining patterns. The latter are due to reactions to several different nuclear antigens.

- Homogeneous:** The entire nucleus fluoresces evenly with a diffuse staining pattern.
- Nuclear Membranous:** The nuclear membrane stains most intensely with decreasing staining intensity of the nucleoplasm towards the center of the nucleus.
- Speckled:** Discrete coarse to fine round speckles fluoresce throughout the nucleus.
- Nucleolar:** The nucleoli stain as multiple solid bodies within the nucleus.
- Centromere:** Large speckles of finite number. Reactive antigen segregates with condensed chromosomes in cells undergoing mitosis.

The specificity of some of the antibodies giving the above staining patterns may be further identified by tests for antibodies to nDNA and to various extractable nuclear antigens. These may be of diagnostic significance as listed in Table 1 at the end of this document.

AMA may be observed on both the distal and proximal tubules of the kidney with the distal tubules staining more brightly. Even though the cytoplasm of the gastric parietal cells also stains, AMA should be quantitated on the kidney.

Staining of the stomach muscularis and kidney glomeruli may also be observed with ASMA, but only ASMA seen on the blood vessel walls of the kidney should be reported.

LKM antibodies exhibit a characteristic granular cytoplasmic staining of the liver parenchymal cells and proximal tubules of the kidney. The reactions on stomach are usually negative.

LIMITATION OF THE PROCEDURE

In some cases, sera positive for ANA may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titers determined. In some cases the presence of two or more antibodies in a serum which are reactive with the same substrate may cause an interference in their detection by immunofluorescence. This interference may cause either failure to detect ANA or suppression of its titer if the interfering antibody has a higher titer than ANA. All ANA reactions should be reported. A positive ANA should not be considered diagnostic of SLE by itself. They also occur in patients with other connective tissue diseases and certain drugs such as procainamide and hydralazine may induce a positive ANA¹. Moreover, sera of patients with malignancies and infectious diseases may also have positive ANA. The clinician should consider the results of all positive indirect immunofluorescence tests along with the results of other laboratory tests and the clinical condition of the patient when making a diagnosis.

EXPECTED VALUES

As seen in Tables 1-7 at the end of this document, tests for nuclear antibodies are used to screen for SLE and certain other immunologic disturbances. AMA occur in over 90% of cases of primary biliary cirrhosis and 3-11% of cases of chronic hepatitis. ASMA occur in the majority of cases of chronic active hepatitis and AGPA are commonly associated with pernicious anemia and chronic atrophic gastritis. Anti-LKM antibodies are found in a subgroup of patients with idiopathic autoimmune chronic hepatitis^{18,20}. The majority of patients with LKM antibodies

EN

had clinical liver disease. Some LKM antibody patients may have sub-clinical hepatitis or hepatocellular carcinoma. LKM antibodies identify a subgroup of HBsAg-negative chronic active hepatitis in which other autoantibody markers are absent. In addition these patients exhibit low levels of serum IgA. See Table 7.

PERFORMANCE CHARACTERISTICS

The ImmuGlo™ Autoantibody Test System was compared with another commercially available fluorescent antibody test using HEp-2 cells as a substrate. The comparison included 15 serum samples from normal subjects as well as sera from patients with the diagnosis of SLE, subacute cutaneous lupus erythematosus, scleroderma or rheumatoid arthritis. Sera were tested according to the procedure and screening dilution recommended by the manufacturer. The results appear in Table 8 at the end of this document.

Sera obtained from 96 normal subjects, 21 patients with SLE, 17 patients with scleroderma and 20 patients with rheumatoid arthritis were tested on the antinuclear antibody (ANA) test (mouse liver sections) kit and other commercially available ANA kits. Sera were tested according to the procedure and screening dilution recommended by the manufacturers. These yielded comparable results as indicated in Table 9.

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REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Tan EM. Autoantibodies to nuclear antigens (ANA): Their immunobiology and medicine. *Adv Immunol* 33: 167-240, 1982.
2. Kumar V, Beutner EH and Chorzelski TP. Autoimmunity of the skin. In "Concepts in Immunopathology", Vol 1, Cruse JM and Lewis RE Jr, Eds, Karger, Basel, 318-353, 1985.
3. Reimer G, Cornell RC and Tan EM. The biochemical nature of nuclear antigens reactive with antinuclear antibodies. In "Immuno- pathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 519-531, 1987.
4. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Standardization of antinuclear antibody and other immunofluorescent tests used in immunopathologic studies of the skin. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 41-64, 1987.
5. Tan EM, Chan EKL, Sullivan KF and Rubin RL. Antinuclear antibodies (ANAs): Diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity. *Clin Immunol Immunopathol* 47: 121-141, 1988.
6. Manns M, Gerken G, Kyriatsoulis A and Meyer zum Büschenfelde KH. Significant autoimmune markers of autoimmune liver disorders: Current status. *J Clin Lab Anal* 1: 362-370, 1987.
7. Mackay IR. Autoimmunity and the liver. *Clin Aspects Immunity* 2: 8- 17, 1988.
8. McMillan SA, Alderdice JM, McKee CM et al. Diversity of autoantibodies in patients with anti-mitochondrial antibody and their diagnostic value. *J Clin Path* 4: 232-236, 1987.
9. Gershwin ME, Coppel RL and Mackay IR. Primary biliary cirrhosis and mitochondrial autoantigens - insights from molecular biology. *Hepatology* 8: 147-151, 1988.
10. Berg PA and Klein R. Mitochondrial antigens and autoantibodies from anti-M1 to anti-M9. *Klin Wochenschr* 64: 897-909, 1986.
11. Popper H and Paronetto F. Clinical, histologic and immunopathologic features of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 339-354, 1980.
12. Berg PA and Bacon H. Serology of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 355-373, 1980.
13. Anderson P, Small JV and Sobieszek A. Studies on the specificity of smooth muscle antibodies. *Clin Exp Immunol* 22: 22-29, 1975.
14. Kurki P, Miettinen A, Linder E, Pikkarainen P, Vuaristo M and Salaspuro MP. Different types of smooth muscle antibodies in chronic active hepatitis and primary biliary cirrhosis: Their diagnostic and prognostic significance. *Gut* 21: 878-884, 1980.
15. Fisher JB and Taylor KB. The significance of gastric antibodies. *Brit J Haematol* 20: 1-7, 1971.
16. Chisholm M. Immunology of gastritis. *Clin Gastroenterol* 5: 419-428, 1976.
17. Bigazzi PE, Burek CL and Rose NR. Antibodies to tissue-specific endocrine, gastrointestinal and neurological antigens. In "Manual of Clinical Laboratory Immunology". Rose NR, Friedman H and Fahey JL, Eds, American Society for Microbi- ology, Washington DC, 762-770, 1986.
18. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopa- thology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
19. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1999 [HHS Pub. No. (CDC) 93-8395].
20. Nisengard RJ. Antinuclear antibodies: Significance of titers. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Bean S, Eds, John Wiley and Sons, New York, 2nd Ed, 387-398, 1979.

21. Meyer zum Büschenfelde KH, Manns M and Trautman F. Autoimmunity in chronic liver diseases - relationship to SLE? In "Recent Advances in Systemic Lupus Erythematosus". Lambert PH, Perrin L, and Izui S, Academic Press, New York, 259-269, 1984.
22. Walker JG, Doniach D, Roitt IM and Sherlock S. Serologic tests in the diagnosis of primary biliary cirrhosis. *Lancet* 1: 827, 1965.
23. Paronetto F and Popper H. Hetero-iso- and autoimmune phenomena in the liver. In "Textbook of Immunopathology", Miescher PA and Müller- Eberhard HJ, Eds, Grune and Stratton, New York, 2nd Ed, 789-817, 1976.
24. Leung PSC, Manns MP, Coppel RL, Gershwin ME. Detection of anti-mitochondrial antibodies in primary biliary cirrhosis and liver-kidney microsomal antibodies in autoimmune hepatitis. In "Manual of Clinical Laboratory Immunology", Rose NR, Hamilton RG, Detrick B, Eds, ASM Press, Washington DC, 6th Ed, 1023-1031, 2002.
25. Muratori P et al. Smooth muscle antibodies and type 1 autoimmune hepatitis. *Autoimmunity*. 35 (8): pp. 497–500. 2002.
26. Gatselis NK et al. Autoantibodies in HCV-treated patients. *World J Gastroenterol*. 11(4):482-487.2005.
27. Miller MH et al. Clinical comparison of cultured human epithelial cells and rat liver as substrates for the fluorescent antinuclear antibody test. *J Rheumatol*. 12 (2): 265-9. 1985.

Table 1. Diagnostic Significance of Antinuclear Antibodies

IF Staining Pattern	Nature of Antigen	Associated Disease
Homogeneous	dsDNA/Histones	SLE
Nuclear Membranous	Laminins	SLE, vasculitis or chronic Hepatitis
Speckled	RNP	SLE or MCTD*
	Sm	SLE
	SS-A/SS-B	SLE or Sjögren's Syndrome
	Scl-70	Scleroderma
Nucleolar	RNAPI	Scleroderma
	Pm-Scl	
Centromere/Kinetochores	inner and outer plates of kinetochores	CREST syndrome

*Mixed Connective Tissue Disease

Table 2: Incidence of Antinuclear Antibodies (ANA) Detected by Indirect Immunofluorescence on HEp-2 Cells

Clinical Condition	No. of Sera	% Positive
SLE	12	100
Subacute Cutaneous LE (SCLE)	7	86
Scleroderma	6	100
Rheumatoid Arthritis	10	50
Normal Controls	15	0

Table 3: Incidence of Antinuclear Antibodies (ANA) Detected by Indirect Immunofluorescence on Mouse Kidney Substrate

Clinical Condition	No. of Sera	% Positive
SLE	21	95
Scleroderma	17	82
Rheumatoid Arthritis	20	5
Normal Controls	96	0

Table 4: Incidence of Anti-Mitochondrial Antibodies (AMA) Detected by IFA on Mouse Kidney Substrate

Clinical Condition	% Positive
Primary Biliary Cirrhosis	100
Autoimmune Chronic Active Hepatitis	8
HBsAg and Chronic Active Hepatitis	0
Extrahepatic Jaundice and Other Liver Diseases	0
Systemic Lupus Erythematosus (SLE)	3
Rheumatoid Arthritis	0
Normal Controls	0

Adapted from Meyer zum Büschenfelde KH, et al.²¹; Walker JG, et al.²² and Paronetto F and Popper H²³.

Table 5: Incidence of Anti-Smooth Muscle Antibodies (ASMA) as Detected by IFA on Mouse Kidney Substrate

Clinical Condition	% Positive
Chronic Active Hepatitis (Type A)	50-87
Primary Biliary Cirrhosis	25
Acute Viral Hepatitis	87
Infectious Mononucleosis	87
Burkitt's Lymphoma	73
Nasopharyngeal Carcinoma	23
Hodgkin's Disease	23
Myeloproliferative Disorder	5
Warts	4
Normal Controls	3-18

Adapted from Anderson P, et al.¹³ Muratori P, et al²⁵

Table 6. Incidence of Anti-Gastric Parietal Cell Antibodies (AGPA) as Detected by Indirect Immunofluorescence on Mouse Stomach Substrate

Clinical Condition	% Positive
Pernicious Anemia (PA)	85-95
Chronic Atrophic Gastritis without PA	30-60
Gastric Ulcer	25-30
Autoimmune Endocrinopathies	25-33
Sjögren's Syndrome	30
First Degree Relatives of PA Patients	30
Normal Controls	
< 20 years old	2
20-60 years old	6-8
> 60 years old	16

Table 7. Characteristics of LKM Antibody Positive Patients

	Condition	% Incidence
Sex (F/M)	4:1	
Age	Mean 34 years	
HBsAg	Negative	
	Other Autoantibodies	
	Anti-nuclear	12%
	Anti-smooth muscle	18%
	Anti-mitochondrial	0%
	Anti-thyroglobulin	12%
	Anti-microsomal	35%
	Immunoglobulin levels (gm/dl)	
	IgG	Normal
	IgM	Normal
	IgA	Low

Adapted from Leung PSC, et al.²⁴

Table 8. Comparison of IFA Antibody Detection Systems

Clinical Condition	No. of Sera	% Positive	
		Immco™	Other
SLE	12	100	100
Subacute Cutaneous LE (SCLE)	7	85	85
Scleroderma	6	100	100
Rheumatoid Arthritis	10	50	30
Normal Controls	15	0	0

Table 9. Comparison of Kits Using Tissue Sections or HEp-2 Cells for the Detection of ANA

Clinical Condition	n	% Positive		
		Immco™ Mouse Liver	Other Mouse Kidney	Other HEp-2
SLE	21	95	95	86
Scleroderma	17	82	88	77
RA	20	25	25	5
Normal Controls	96	0	0	0

Adapted from Miller MH, et al.²⁷

For technical assistance please contact:

IMMCO Diagnostics, Inc.

60 Pineview Drive

Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immco.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

www.emergogroup.com