Anti- β_2 GP1 Antibody ELISA

IVD

CLIA Complexity: High
 CDC Analyte Identification Code: 0529
 CDC Test System Identification Code: 28484 (IgG), 28494 (IgM)

PRODUCT INSERT

REF 1152G Anti- β_2 -GP1 IgG ELISA 96 Determinations

REF 1152M Anti- β_2 -GP1 IgM ELISA 96 Determinations

An enzyme linked immunoassay (ELISA) for the detection and semi-quantitation of IgG and IgM antibodies to β_2 -GP1, as an aid in assessing the risk of thrombosis in patients with Systemic Lupus Erythematosus (SLE) or lupus like disorders.

SUMMARY AND EXPLANATION

Antiphospholipid antibodies are a heterogeneous group of autoantibodies against negatively charged phospholipids¹. They are detected primarily by the anti-cardiolipin antibody (ACA) test, the biological false positive test for syphilis or the lupus anticoagulant test. These three tests detect related, but not necessarily identical antibodies. Thus, more than one of these tests may be necessary to identify antiphospholipid antibodies.

Of the various tests for the detection of antiphospholipid antibodies, the anti-cardiolipin antibody test performed by ELISA is the most sensitive². The presence of anti-cardiolipin antibodies helps to identify patients at risk of venous and/or arterial thrombosis often accompanied by *thrombocytopenia*, a syndrome referred to as *antiphospholipid syndrome*¹⁻¹². It most commonly occurs in patients with *systemic lupus erythematosus* (SLE) or lupus-like disease where the criteria for SLE are not fulfilled⁵⁻⁷. High levels of anti-cardiolipin antibodies occur in thrombosis, fetal loss, *thrombocytopenia* and several other disorders¹⁻¹⁵. Anti-cardiolipin antibodies are routinely detected by ELISA, using cardiolipin as the antigen. Recent observations have indicated that 50kD serum factor is necessary for ACA to bind to cardiolipin. This co-factor was later identified as β_2 -GP1 (β_2 -glycoprotein I) or synonymously apolipoprotein H¹⁶⁻¹⁸. Though the function of β_2 -GP1 remains unclear it is certain that the presence of β_2 -GP1 facilitates the binding of the ACA to cardiolipin. Anti-cardiolipin antibodies, especially at low levels are found in a variety of clinical disorders unrelated to anti-phospholipid syndrome. ACA present in sera of patients with SLE and other autoimmune disorders is immunologically distinct from sera of patients with syphilis. The use of β_2 -GP1 protein as a replacement for cardiolipin antigen on the solid matrix, helps to make this distinction. Furthermore, anti- β_2 -GP1 antibodies are very specific for the anti-phospholipid syndrome.

PRINCIPALS OF THE PROCEDURE

The test is performed as a solid phase immunoassay (ELISA) in β_2 -GP1 coated microwells. Controls, Calibrators and patient serum samples are incubated in the antigen coated microwells to allow antibodies present in the serum to bind. Unbound antibody and other serum proteins are removed by washing the microwells. Antibodies bound to the microwells are detected by adding an enzyme labeled anti-human IgG or IgM conjugate to the microwells. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the appropriate class. Unbound enzyme-labeled conjugate is removed by washing. Specific enzyme substrate (pNPP) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of the pNPP substrate. The reaction is stopped and the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in enzyme units per milliliter (EU/ml).

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials¹⁹. **WARNING** - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this kit insert to en-sure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc.

Use good laboratory techniques to minimize microbial and chemical contamination. Do not use after expiration date.

Materials provided

Immulin™ β₂-GP1 IgG ELISA [REF] 1152G

Immulin™ β₂-GP1 IgM ELISA [REF] 1152M

Kits contain sufficient reagents to perform 96 determinations each.

12 x 8	MICROPLATE B2GP1	Microplate with individual breakaway microwells coated with β ₂ -GP1 antigen.
1 x 1.5 ml	CALIBRATOR A B2GP1 *†	Ready to use Calibrator A (<i>green cap</i>). Human serum containing antibodies to β ₂ -GP1.
1 x 1.5 ml	CALIBRATOR B B2GP1 *†	Ready to use Calibrator B (<i>violet cap</i>). Human serum containing antibodies to β ₂ -GP1.
1 x 1.5 ml	CALIBRATOR C B2GP1 *†	Ready to use Calibrator C (<i>blue cap</i>). Human serum containing antibodies to β ₂ -GP1.
1 x 1.5 ml	CALIBRATOR D B2GP1 *†	Ready to use Calibrator D (<i>yellow cap</i>). Human serum containing antibodies to β ₂ -GP1.
1 x 1.5 ml	CONTROL + B2GP1 *†	Ready to use Positive Control (<i>red cap</i>). Contains human serum positive for β ₂ -GP1.
1 x 1.5 ml	CONTROL - *	Ready to use Negative Control (<i>whitecap</i>). Contains human serum.
1 x 12 ml	IgG-CONJ ALKPHOS *†	Ready to use anti-human Alk. Phos. Conjugate . Color coded pink.
1 x 12 ml	IgM-CONJ ALKPHOS *†	Ready to use anti-human Alk. Phos. Conjugate . Color coded pink.
1 x 60 ml	DIL *	Ready to use Serum Diluent . Color coded blue.
1 x 12 ml	SUBSTRATE *	Ready to use Enzyme Substrate . Contains pNPP. Protect from light.
1 x 12 ml	STOP	Ready to use Stop Solution .
2 x	BUF WASH	Powder Wash Buffer . Reconstitute to one liter each.

* Contains <0.1% NaN₃

EN

† **REF** 1152G contains **B₂-GP1** IgG calibrators, controls and IgG conjugate

REF 1152M contains **B₂-GP1** IgM calibrators, controls and IgM conjugate

Symbols used on labels:

LOT Lot number

REF Catalog number

 Use by

 Storage temperature

 Read instructions for use

IVD In vitro diagnostic use

 Manufacturer

 Number of Tests

Materials Required But Not Provided

Pipette capable of dispensing 5 μ l and 1000 μ l

- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Deionized or distilled water
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Squeeze bottle to hold diluted wash buffer
- Timer
- Absorbent paper

Automatic microplate washer capable of dispensing 200 μ l

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.


PROCEDURE

Procedural Notes

- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-26°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

Assay procedure

- ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-26°C) PRIOR TO BEGINNING THE ASSAY.**
- Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.
- Qualitative determination:** use only Calibrator D.
Semi-quantitative determination: use Calibrators A - D as shown in the ex-ample below.

- Prepare a 1:101** dilution of the patient specimen by mixing 5 μ l of the patient specimen with **0.5 ml** of Serum Diluent.
- Add **100 μ l** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record.
Note: Include one well with **100 μ l** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of this well should not be greater than 0.3.
- Incubate **30 minutes** (\pm 5 minutes) at room temperature on a level surface.
- Wash step: Thoroughly aspirate the contents of each well. Add 200-300 μ l of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
- Add 100 μ l of the Conjugate to each well.
- Incubate the wells for **30 minutes** (\pm 5 min) at room temperature.
- Wash step: Repeat step 7.
- Add 100 μ l of Enzyme Substrate to each well.
- Incubate the wells for **30 minutes** (\pm 5 min) at room temperature.
- Add 100 μ l of Stop Solution to each well. Maintain the same sequence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
- Read the absorbance (OD) of each well at 405nm using a single or dual wave-length (405/630nm) microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3 . Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The Negative Control must be <20 EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml. When

EN

performing qualitative determinations, the optical density of Calibrator D must be greater than that of the Negative Control and less than the absorbance of the Positive Control. For semi-quantitative determinations, the Positive Control must give values in the range stated on the vial.

RESULTS

Calculations

The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

2. SEMI-QUANTITATIVE DETERMINATION

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in EU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value.

Calibrator

The calibrators are included to provide semi-quantitation and must be used with each test. Patient specimens containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml, multiply the units obtained by the dilution factor. See sample standard curves (Figure 1) at the end of this document.

Interpretation

The information in Figure 2 at the end of this document serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

LIMITATIONS OF THE PROCEDURE

The ImmuLisa™ Anti - β_2 -GP1 test should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only.

Furthermore, a diagnosis cannot be made on the basis of anti - β_2 -GP1 results alone. The results of other laboratory tests and clinical findings must be considered. When a negative anti - β_2 -GP1 test occurs in the presence of clinical indications, a lupus anticoagulant test or other additional testing is recommended.

Expected Values

The incidence of anti β_2 -GP1 in SLE with thrombosis and pregnancy is summarized in the following table:

Incidence of ACA and Anti β_2 -GP1 Antibodies in SLE¹⁷

	ACA		Anti β_2 -GP1	
	IgG %	IgM %	IgG %	IgM %
SLE	10	6	12	4
with thrombosis	24	14	27	14
no thrombosis	8	5	10	3
in pregnancy	12	7	15	7
with recurrent fetal loss	20	0	27	0
no fetal loss	11	7	14	7

PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ anti- β_2 -GP1 IgG and IgM Antibody ELISA tests was determined by comparing the results with other commercially available anti- β_2 -GP1 IgG and IgM ELISA methods.

Normal Range: The normal range was established by testing 64 serum samples from apparently healthy donors obtained from the Red Cross. The mean plus three standard deviations of the mean of this normal population was used to determine the cut-off between normal and borderline positive individuals.

Comparative Specificity and Sensitivity:

Test results for the comparison between the Immulisa™ anti- β_2 -GP1 IgG and IgM Antibody tests and other commercially available assays for the same analytes appear in Tables 1 and 2 at the end of this document.

Cross-reactivity:

It is known that patients suffering from infectious diseases like syphilis have cardiolipin antibodies and lupus anti-coagulant. To test the clinical specificity of the assays, we tested sera from patients with syphilis for the presence of anti-cardiolipin and anti- β_2 -GP1 antibodies. The results appear below.

	# Tested	# Positive	% Positive
Anti-Cardiolipin (ACA) IgG	49	49	100
Anti- β_2 -GP1 IgG	49	2	4
Anti-Cardiolipin (ACA) IgM	50	17	34
Anti- β_2 -GP1 IgM	50	4	8

These studies suggest that β_2 -GP1 antibodies are more specific than anti-cardiolipin antibodies as only a small percentage of the patients with syphilis are positive for β_2 -GP1 antibodies.

Precision:

Two samples with known concentrations of anti β_2 -GP1 were assayed in 10 replicates over a two weeks period. Intra-and inter-assay coefficient of variation (CV) appear in tables 3 and 4 at the end of this document.

Recovery:

Three samples with known concentrations of anti - β_2 -GP1 were mixed with the appropriate dilutions of another positive sample with known levels of anti - β_2 -GP1. The anti - β_2 -GP1 values of the spiked samples were then determined and from the values obtained the percent recovery calculated. The results appear in tables 5 and 6 at the end of this document.

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REFERENCES I ΛΙΟΓΡΑΦΙΑ LITERATUR BIBLIOGRAPHIE BIBLIOGRAFIA

1. Harris EN, Phil M, Asherson RA et al. Antiphospholipid antibodies - antibodies with a difference. *Ann Rev Med*; 1988, 39:261-271.
2. Harris EN, Gharavi AE, Wasley GD and Hughes GRV. Use of an enzyme-linked immunosorbent assay and inhibition studies to distinguish between antibodies to cardiolipin from patients with syphilis or autoimmune disorders. *J Infect Dis*; 1988, 157:23-31.
3. Lockshin MD. Anti-cardiolipin antibody. *Arth Rheum*; 1987, 30:471-472.
4. Derksen RHW, Beisma D, Bouma BN et al. Discordant effects of prednisone on anticardiolipin antibodies and the lupus anti-coagulant. *Arth Rheum*; 1986, 29:1295-1296.
5. Alarconegovia D, Deleze M, Oria CV et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. *Medicine*; 1989, 68:353-365.
6. Triplett DA. Antiphospholipid antibodies and recurrent pregnancy loss. *Am J Reprod Immunol*; 1989, 20:52-67.
7. Montecucco C, Longhi M, Caporali R et al. Hematological abnormalities associated with anticardiolipin antibodies. *Hematologica*; 1989, 74:195-204.
8. Asherson RA, Khamashta MA, ORDI-Ros J, et al. The "primary" antiphospholipid syndrome: major clinical and serological features. *Medicine*; 1989, 68:366-374.
9. Walport MJ. Pregnancy and antibodies to phospholipids. *Ann Rheum Dis*; 1989, 48:795-797.
10. Machworth-Young CG, Loizou S and Walport MJ. Antiphospholipid antibodies and disease. *Medicine*; 1989, 72:767-777.
11. Harris EN, Gharavi AE, Hughes GRV. Antiphospholipid antibodies. *Clin Rheumatol*; 1985, 1:591-609.
12. Lubbe WF, Palmer SJ, Butler WS and Laggins GC. Fetal survival after prednisone suppression of maternal lupus anticoagulant. *Lancet*; 1983, 1:1361-1363.
13. Hogan MJ, Brunet DG, Ford PM et al. Lupus anticoagulant, antiphospholipid antibodies and migraine. *Can J Neurological Sci*; 1988, 15:420-425.
14. Harris EN, Pierangeli S. Antiphospholipid antibodies: method of detection. *Am J Reproduct Immunol*; 1992, 28:208-210.
15. Silveira LH, Lopez LR, Uzzle RT et al. IgA, IgG and IgM anticardiolipin antibodies in black patients with systemic lupus erythematosus. *Arthr Rheum 56th Annual Meeting*; 1992, S125.
16. McNeil HP, Simpson RJ, Chesterman CN and Krilis SA. Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation: β_2 -glycoprotein I (apolipoprotein H), *Proc Natl Acad Sci USA*. 87:4120-4124, 1990.
17. Matsuura E, Igarashi Y, Jugimoto M, Ichikawa K, Suzuki T, Sumida T, Yasuda T and Koike T. Heterogeneity of anticardiolipin cofactor. *J Immunol*. 148:3885-3891, 1992.
18. Matsuura E, Igarashi M, Igarashi Y, Nagae H, Ichikawa K, Yasuda T and Koike T. Molecular definition of human β_2 -glycoprotein I (β_2 -GPI) by cDNA cloning and inter-species differences of β_2 -GPI in alternation of anticardiolipin binding. *Int Immunol*. 3:1217-1221, 1991.
19. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health; 1993, (HHS Pub. No [CDC] 93-8395).

Figure 1. Standard Curve


Figure 2. Interpretation

Anti- β_2 -GP1 IgG/IgM	
Value	Interpretation
<20 EU/ml	Neg (-)
20 - 25 EU/ml	Borderline
>25 EU/ml	Pos (+)

Table 1. ImmcoB Anti- β_2 -GP1 IgG vs. Other ELISA

		ImmcoB AGA-IgA		
		POS (+)	NEG (-)	TOT (=)
Other ELISA	POS (+)	17	0	17
	NEG (-)	1	74	75
	TOT (=)	18	74	92

relative specificity: 99%
 relative sensitivity: 100%
 relative agreement: 99%

Table 2. ImmcoB Anti- β_2 -GP1 IgM vs. Other ELISA

		ImmcoB AGA-IgG		
		POS (+)	NEG (-)	TOT (=)
Other ELISA	POS (+)	16	2	18
	NEG (-)	1	49	50
	TOT (=)	17	51	68

relative specificity: 96%
 relative sensitivity: 89%
 relative agreement: 98%

Table 3. Anti- β_2 -GP1 IgG Precision

β_2 -GP1 IgG	Intra-plate	Inter-plate
	CV	CV
1.	7.8%	6.2%
2.	9.7%	7.0%

Table 4. Anti- β_2 -GP1 IgM Precision

β_2 -GP1 IgM	Intra-plate	Inter-plate
	CV	CV
1.	8.1%	11.4%
2.	9.6%	11.7%

Table 5. Anti- β_2 -GP 1 IgG Recovery

IgG	EU/ml added	EU/ml obtained	% Recovery
1.	109.5	114.8	105%
2.	55	59	107%
3.	45.5	45.1	99%

Table 6. Anti- β_2 -GP 1 IgM Recovery

IgM	EU/ml added	EU/ml obtained	% Recovery
1.	169.5	168.7	100%
2.	84.5	87.2	103%
3.	70.5	75.7	107%

For technical assistance please contact:

IMMCO Diagnostics, Inc.

60 Pineview Drive

Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immco.com

or your local product distributor

EU A thori ed Representati e/A torisierter Repr sentant/Rappresentante

A torri ato/Representante A tori ado/Repr sentant A toris

EMERGO Gro p, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

www.emergogroup.com